

WinnerPage 蛋白预制胶 the 🏆inner, the better!

■ 产品简介

WinnerPage 蛋白预制胶是一款安全、快捷、高性能的预制聚丙烯酰胺凝胶，可用于蛋白质分离。本产品提供不同浓度的梯度胶和固定浓度胶，并有12孔和15孔两种孔数选择。采用全自动凝胶灌注技术，产品的重复性好，质量稳定。独特的凝胶缓冲配方使蛋白电泳条带更为清晰锐利，更加均匀，分辨率更高。配套缓冲液为中性缓冲液，可以提高凝胶稳定性和避免蛋白在电泳过程中的再修饰。

■ 产品特点



操作方便

即开即用，试剂盒中**提供配套缓冲液**，
无需自己配制各种缓冲液和灌胶操作；



分辨率高

全新凝胶缓冲体系配方使蛋白电泳条带
更清晰锐利，分辨率更高；



节省时间

最快30 min完成电泳；



兼容性强

兼容市场上主流的mini电泳槽，如
Bio-Rad、天能、六一等；



安全无毒


无需接触有毒试剂；



稳定性好

大规模自动化生产技术，保证产品良
好的重复性和质量的稳定性；

不同浓度预制胶分离图谱 (Tris/MOPS/SDS 电泳缓冲液, 蛋白条带分子量单位: kDa)



产品信息

目录号	产品名称	规格	目录价(元)
CW8201S	WinnerPage Protein Gel, 8%, 12 wells, Bis-Tris	10 PCs/Box	268
CW8202S	WinnerPage Protein Gel, 8%, 15 wells, Bis-Tris	10 PCs/Box	268
CW8203S	WinnerPage Protein Gel, 10%, 12 wells, Bis-Tris	10 PCs/Box	268
CW8204S	WinnerPage Protein Gel, 10%, 15 wells, Bis-Tris	10 PCs/Box	268
CW8205S	WinnerPage Protein Gel, 12%, 12 wells, Bis-Tris	10 PCs/Box	268
CW8206S	WinnerPage Protein Gel, 12%, 15 wells, Bis-Tris	10 PCs/Box	268
CW8207S	WinnerPage Protein Gel, 4-12%, 12 wells, Bis-Tris	10 PCs/Box	268
CW8208S	WinnerPage Protein Gel, 4-12%, 15 wells, Bis-Tris	10 PCs/Box	268
CW8209S	WinnerPage Protein Gel, 4-20%, 12 wells, Bis-Tris	10 PCs/Box	268
CW8210S	WinnerPage Protein Gel, 4-20%, 15 wells, Bis-Tris	10 PCs/Box	268
CW8211S	MOPS-SDS Running Buffer Powder	10 PKs/Box	208



4006-222-360 (免费电话)
service@cwbiotech.com
www.cwbio.com



让每一个生命健康有为
上海证券交易所上市公司 (688426.SH)